HNK-SUTPUMP-01

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Index

1.	Introduction.....	23
2.	Device Specifications	23
3.	Technical Specifications.....	24
4.	Design and Specifications	24
5.	Important Information.....	26
6.	Information for Use.....	28
7.	Preparing the Product for the Use	29
8.	Breast Pump Use	31
9.	After Use	33
10.	Power Options	34
11.	Maintenance	35
12.	Storage and Transportation	35
13.	Troubleshooting.....	36
14.	Warranty Conditions and After Sale Support.....	37
15.	EMC (Elektromagnetic Compatibility) Decleration	37
16.	Symbols and Descriptions.....	41
17.	Device & Outside Box's Symbols Descriptions.....	42
18.	Manufacturer and Technical Service Information	42

1. Introduction

Thank you for choosing **Nimo**® Electronic Single Breast Pump. HNK-SUTPUMP-01 is designed to be ergonomic and easy to use. It has an inclined funnel design that will facilitate suction while sitting in an upright position and a chest pad to give your chest comfort.

It also has a silicone diaphragm to prevent milk from escaping to the tube and back to the pump. The funnel, chest cushion, silicone diaphragm, silicone valve and baby bottle can be washed in the dishwasher.

The breast pump starts working in low-level mode, helping to start the mother's milk flow. The mother can choose the optimal rhythm of milking with the 4-stage suction level. The milk pump has two different suction modes, which can be absorbed by the bionic rhythm. The device has also a classical fixed rhythms that mimic the suction rhythms of a baby. Both modes have 4 different suction power levels. Thanks to the LED indicators, it is easily seen in which suction mode and level the device is operating.

It can be used with 4 AA batteries, in addition it is possible to use the device with electrical connection.

2. Device Specifications

- Milking with battery or adapter.
- A bionic rhythm, mimicking the infant sucking
- Soft silicone cushion for painless and comfortable use
- The device has 8 different rhythm options, also 4 different suction power levels which is adjustable in fixed rhythm and bionic rhythm modes.
- Funnel design that allows manual operation
- Quiet and vibration-free operation
- Automatically shut down after 30 min.
- 150 ml bottle capacity
- Easy to use and ergonomic design


3. Technical Specifications

Product Name	Nimo ® Electronic Single Breast Pump
Model Number	HNK-SUTPUMP-01
Max. Vacuum	160 mmHg
Operating Conditions	5 ~ 40 ° C ≤ 95 % RH
Storage & Delivery Conditions	Temperature: -25 ~ 55° C Humidity: ≤ 95 % RH Atmospheric pressure: 70 ~ 106 kPa
Power	6W-DC 6V 1A veya (4x1.5 AA Alkaline Batteries) (Batteries are not included !)
Adapter	Input: 100-240VAC ~ 50/60Hz Output: DC 6V 1A
Automatically Shut Down Time	30 min
Solid and Liquid Protection Class	IP22
Classification	Class II
Sizes	105mm(length) x 100mm(width) x 60mm(height)
Weight	225 gr (without Batteries)

4. Design and Specifications

a) Device Image and Parts


- 1- Control and motor unit
- 2- AC/DC Adapter
- 3- Bottle
- 4- One-way milk valve
- 5- Silicone Breast Pad
- 6- Silicone suction diaphragm
- 7- Pump body
- 8- Funnel cover
- 9- Air tube
- 10- Air tube-motor unit connection

b) Button Specifications

Button	Function
	Switching device on and off
	Change the milking mode
	Increase the vacuum level
	Decrease the vacuum level

c) Led Level-Mode Indicator and Meanings

5. Important Information

Please read and keep the user manual carefully before use.

For your child's safety and health WARNING!

- ! Always use this product with adult supervision.
- ! Never use the feeding teats as a soother.
- ! Continuous and prolonged sucking of fluids will cause tooth decay.
- ! Always check food temperature before feeding.
- ! Check the bottle, teat and all parts by pulling them in all directions before each use. Do not use if you see signs of damage or thinning.
- ! Wash your hands thoroughly and make sure that the surfaces of the sterilized parts which will get in body contact are clean.
- ! Before first use place in boiling water for 5 minutes.
- ! At the first use and later after each use, remove all parts of the device and sterilize it by using a sterilizer for hygiene or by boiling in hot water for 5 minutes.
- ! Do not expose teat to the disinfectant, direct sunlight or direct heat for longer than recommended, as this may cause wear.
- ! Do not place in a heated oven.
- ! Dishwasher safe, food colorings may discolour components and scratches may occur on the product.
- ! DO NOT use abrasive cleaning detergents or anti-bacterial cleaning products. The use of concentrated detergent may cause deterioration or cracking of the bottle parts. Avoid this and do not use the product if this happens.
- ! DO NOT heat the baby food and the bottle in the microwave. The food in the baby bottle may not heat equally, it can cause your baby's mouth to burn.

- ❗ Because of hygiene, you should not increase the holes of the bottle pacifiers yourself. As your baby grows, use the appropriate flowing baby bottle. It is recommended that you replace the bottle pacifier with a new one after 1-2 months of use.
- ❗ It is not recommended to use fruit juices and flavored sugary drinks except than milk and water.
- ❗ Keep unused parts out of the reach of children. DO NOT allow your child to play with small parts or to walk / run while using the bottle.

ATTENTION

- If the pump is not working properly, do not use it and contact the supplier.
- Do not use the pump during pregnancy.
- The appliance is hygienic, so use the device only for yourself.
- If you feel pain or discomfort during use of the pump, stop the use and consult your doctor.

WARNINGS

- Do not expose the motor unit and its electronics to water or moisture, clean with a wet cloth, then wipe with a dry cloth.
- Operate the **Nimo**® Electronic Single Breast Pump with the original adapter or 4x1.5V AA Alkaline battery. After each use, check the device is switched off and unplug it from the wall outlet.

- Maximum milk level should not be higher than the bottle capacity (150 ml).
 - Always hold the bottle upright.
-

TIPS

- For mothers who will use the device for the first time, it is recommended to stimulate the breast tips for 2-5 minutes before use.
- Massaging your breasts before milking allows you to get better milk.
- Your breast tips should not be in contact with the breast pump head.

6. Information for Use

BREAST PUMP USE TIME

Use the breast pump when the milk comes or if you want to store the milk.

Start using the breastpump 2 or 4 weeks after giving birth, unless your doctor recommends otherwise.

OTHER CONDITIONS

- You can use it when you can not feed the baby directly with your breast.
- If you want to feed your baby later, you can store the milk with the help of breast pump.
- If you are separated from your baby and want to continue breastfeeding when you're back together, you must provide your milk regularly to stimulate your milk production.

7. Preparing the Product for the Use

- 1) Install the unidirectional milk valve in the location below the pump housing. Make sure that it is seated in place by pushing it as far as possible.

- 2) Fit the pump body by turning it clockwise until the bottle is seated well.

- 3) Install the suction level adjuster part in the silicone suction diaphragm. Make sure it is seated properly.

- 4) Install the combined silicone suction diaphragm and suction level adjuster on the upper side of the pump body. Make sure that the frame is fully engaged by pressing it down.

- 5) Fit the funnel cover into the pump housing as shown in the picture and make sure that it fits properly.

- 6) Attach the end of the air tube to the inlet of the funnel cover as shown at the picture and make sure that no air escapes.

- 7) Connect the tube and the motor unit. Please connect the end of the air tube as shown at the picture.

- 8) Insert the air tube, the end of which is attached to the clamp, into the air inlet of the motor unit as shown at the picture.

8. Breast Pump Use

- Wash your hands thoroughly. Make sure your breasts are clean.
- Check that the pump is installed correctly.
- Sit in a comfortable upright position. You can use the pillow to support your back. Please do not use your pump while lying down.
- Push the assembled breast pump into your breast. Adjust the breast tip so that the nipple is right in the middle.

- Press  button to turn on the device and to start the vacuum process.

Natural (bionic) Rhythm Mode

- The device starts at level 1 of the **natural (bionic) rhythm mode**. The natural (bionic) rhythm mode has 4 different vacuum levels and mimicks the baby's sucking reflex. In this rhythm, the device massages the mother's breast and makes it easy for the mother to feel comfortable.
- ❗ **Important:** Mothers with breast tip problems or with sensitivities, the device can be used in natural (bionic) rhythm mode.
- Set the desired vacuum level by using   buttons.
- Press  button to switch from the natural (bionic) rhythm to the fixed rhythm mode.
- To switch off the device, press  button.

Fixed Pattern Mode

- The constant rhythm mode of the device is the normal operating rhythm mode. The mode has 4 different vacuum levels and it is recommended to select the best suited vacuum level for the mother.
- In this mode the appropriate vacuum level is determined by using the   buttons
- Press  button to switch from the natural (bionic) rhythm to the steady rhythm mode.
- To turn off the device press  button.

- ! **Important:** The vacuum levels of natural (bionic) and fixed rhythm modes, can be easily monitored from LED indicators.
- ! **Important:** In order to use Nimo® Electronic Single Breast Pumps efficiently, we recommend you to use it first in (bionic) rhythm mode. After your milk has come, we recommend the user to change it in the steady rhythm.
- ! **Warning:** If you are unable to milk your milk, do not pump it for more than 5 minutes. Try to relax and do the milking another time during the day.

9. After Use

STORAGE OF MILK

Your baby milk should be stored in a sterile storage bottle only.

You can keep the milk which comes from your breast at 40 °C in the bottle heater and feed your baby in the first 4 hours.

To freeze the milk, only fill up ¾ of the baby bottle, allow the expansion of the milk
Note the date of the milking at the cover of the bottle and use the old milk first.

Storage Place	Storage Place Temperature	Breast Milk Storage Period
Room Temperature	19-25 °C	6 hours
Refrigerator	>4 °C	3 days
Refrigerator	<4 °C	5 days
Refrigerator Freezer Compartment		14 days
Deep Freezer	<-18 °C	6 months

The milk which you remove from the freezer should be thawed naturally. In case of emergency, the heater can be thawed with the help of the heater.

Thawed milk can be stored in the refrigerator for 24 hours.

- ! **Warning:** Never freeze thawed milk.

Frozen milk can be heated slowly in warm water or in a bottle heater.


- ! **Warning:** Never heat the milk in the microwave. May cause antibodies or food poisoning.
- ! **In addition, milk can cause overheating and boiling. May cause irritation to the throat.**
- ! **Do not place the stored milk, that you removed from the freezer directly in boiling water. It may cause fractures and cracks in the milk container.**
- ! **Do not place boiled water directly into the bottle, let it cool for 20 minutes.**

10. Power Options

AC/DC Adapter (Included in the package)

Connect the AC Adapter to the power connector socket on the right side of the motor unit. Plug the AC Adapter into an empty outlet. Electrical safety is ensured when the breast pump is only plugged into the supplied AC/DC power adapter.

4X1.5V AA Alkaline Battery (Not Included)

Open the battery protection cover located on the backside of the device and gently press down on the battery cover after placing 4 AA alkaline batteries as shown with the + and - terminals in the device's battery bay.

! **Warning:** It is not recommended to use a rechargeable battery.

! **Note:** If you do not use the device for a long time, please remove the batteries from the device. Replace all batteries at the same time.

! **Note:** Dispose of used batteries in accordance with environmental protection and recycling regulations.

11. Maintenance

- ⚠ **WARNING:** After sterilizing the parts of the device, make sure that your hands are clean before contacting with your hands. Make sure that the surfaces in which the parts come into contact during cleaning are clean.
- ⚠ **WARNING:** Remove, clean and sterilize all parts before using the product for the first time.

1. After each use, remove all parts of the pump and wash with hot soapy water and rinse thoroughly. All Parts of the appliance can also be washed in the dishwasher.
2. Clean the motor unit **ONLY** with a wet cloth and then dry it afterwards. Do not use polishing, abrasive or anti-bacterial materials to clean the product.
3. Sterilize the cleaned parts by a sterilizer or boiled in water for 5 minutes.
4. Check that the vacuum tube connections are fully inserted and that the suction valve is not damaged.
5. The pump housing and milk storage bottles can be safely washed in the dishwasher, but may be discolored due to the colors from the food.
6. After rinsing the sterilized parts on a clean surface, install all pieces together and keep them out of reach of children.

NEVER sterilize and clean the air tube connection. Remaining Water in the tube or connection part may damage electronic parts.

12. Storage and Transportation

The device can be stored in a well ventilated room with the package at a temperature of -25 ~ 55 ° C, which will not exceed 95% humidity.

Protect the package from wetting, shocks and crushing during transportation. Position according to the marked direction.

➤ Package Content

- 1 piece HNK-SUTPUMP-01 electronic single breast pump.
- User Manual
- AC/DC Adapter
- Replacement unidirectional suction valve

- Warranty certificate
- Accessories
 - Silicone baby bottle pacifier

13. Troubleshooting

Problem	Solution
Device is not working.	☞ Make sure the air tube is installed in the motor unit or pump housing. ☞ If the device is powered by a battery, make sure that the batteries are correctly inserted and charged. ☞ If the device is powered by an AC / DC adapter, make sure that the adapter is secure and is correctly connected to the device.
The pump vacuum power is too weak or the pump does not vacuum at all.	☞ Make sure that the unidirectional suction valve is not damaged. ☞ Make sure that all connections are correct. Make sure the vacuum diaphragm is seated correctly. ☞ Make sure the breast pump cap is properly seated in your breast. ☞ Push the chest funnel into your chest first, then try to operate the device.
Milk is not reaching the bottle	☞ Make sure your breast tip is centered on the milk pump cap. ☞ Make sure that the unidirectional suction valve is not damaged. ☞ Make sure the milk is coming or not from your breasts

14. Warranty Conditions and After Sale Support

The warranty period of the device is 2 years from the date of purchase.

- ! **Warning:** If the device is opened, damaged or repaired by the user, it is not covered by the warranty.
- ! **Note:** Please do not forget to take and store your sales receipt with warranty certificate.

15. EMC (Elektromagnetic Compatibility) Declaration

This equipment generates, uses and can radiate radio frequency (RF) energy. If this equipment is not installed and used in accordance with the instructions, it may cause electromagnetic interference.

This device has been tested for “Electromagnetic Compatibility – Requirements for Household Appliances, Electric Tools and Similar Apparatus” in accordance with EN-55014-1 and has been determined to comply with acceptable limits. These limits indicate that the device provides acceptable levels of protection against electromagnetic interference (EMC) if the device is used as specified in the manual.

This device may be affected by portable and mobile RF communication devices. This equipment must not be stored with other equipments.

For more information about this device and the EMC, see the tables below

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment – guidance
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	EN 55014-1 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge Immunity Test(ESD) IEC 61000-4-2	± 2 kV, ± 4 kV contact ±2 kV, ±4 kV, ±8 kV air	± 2 kV, ± 4 kV contact ±2 kV, ±4 kV, ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst Immunity Test IEC 61000-4-4	1 kV At 5 kHz, 2 min coupling time	1 kV At 5 kHz, 2 min coupling time	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 1 kV	± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 2 kV	± 2 kV	
Voltage Disp, Short Interruptions and Voltage Variations Immunity Test (Voltage Dips) IEC 61000-4-11	0% <i>UT</i> ;0.5 cycle 10ms 40% <i>UT</i> ;10 cycle 200ms 70% <i>UT</i> ;25 cycle 500ms	0% <i>UT</i> ;0.5 cycle 10ms 40% <i>UT</i> ;10 cycle 200ms 70% <i>UT</i> ;25 cycle 500ms	Mains power quality should be that of a typical commercial or hospital environment.

Not: *UT* is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	EN 55014-1 Test level	Compliance level	Electromagnetic environment – guidance
Immunity to Conducted Disturbances, Induced by Radio Frequency Fields IEC 61000-4-6	150 kHz to 230 MHz 3Vrms 80% Am at 1kHz Sinusoidal 1% step width 2 s dwell time	150 kHz to 230 MHz 3Vrms 80% Am at 1kHz Sinusoidal 1% step width 2 s dwell time	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distance: $d = 0.35\sqrt{P}$ $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ $d = 2.3\sqrt{P}$ 80 MHz to 800 MHz 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: (((•))) ▲
Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands)	150 kHz to 80 MHz (out ISM and amateur radio bands)	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1.2\sqrt{P}$	$d = 0.6\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.06	0.12	0.23
0.1	0.38	0.19	0.38	0.73
1	1.2	0.6	1.2	2.3
10	3.8	1.9	3.8	7.3
100	12	6	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	EN 55014-1 Test Level	Compliance level	Electromagnetic environment – guidance
Discontinuous Disturbance / Click Test	150 kHz, 500 kHz, 1.4 MHz, 30 MHz	150 kHz, 500 kHz, 1.4 MHz, 30 MHz	Device is supplied by LISN equipment on L, N, PE couplings in the shielded room. When device is in the normal period of operation click value that was given to main supply is measured by EMI Test Receiver.
Conducted Emission	150 kHz – 30 MHz, Quasi-Peak, Average Detector.	150 kHz – 30 MHz,	
Power Disturbance Test	30 MHz – 300 MHz	30 MHz – 300 MHz	Device, has been placed on a wooden isolator which was 80 cm height from reference ground plane in the shielded room. It is tested to find the voltage value to produce the maximum emission on the frequencies 30 MHz and 300 MHz supply voltage.

16. Symbols and Descriptions

	CE Marking
	Manufacturer
	Date of Manufacture
	Consult instructions for use
	Suitable for domestic use
	Direct current
IP22	Solid and Liquid Protection Class
LOT	LOT Number
REF	Model Number
SN	Serial Number
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.

17. Device & Outside Box's Symbols Descriptions

	Warning/Caution/Note
	Storage and Handling Ambient Temperature Range: -25 °C ~ 55 °C
	Storage and Handling Environment Humidity Range : ≤95% RH
	Storage and Handling Environment Pressure Range : 70 ~ 106 kPa
	Keep Away From Sunlight
	Keep Away From Rain
	Food Safe Symbol
	Barcode Number

18. Manufacturer and Technical Service Information

Hünkar Ecza ve Medikal San. Tic. Ltd. Şti.
Akçeşme Mahallesi 2019 Sk. No:5 Merkezefendi
Denizli /TURKEY
Phone: +90 258 371 46 56
www.hunkarmedikal.com