HNK-MESH-01

Tařınabilir Mesh Nebulizatör
Portable Mesh Nebulizer

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Hünkar Ecza ve Medikal San. Tic. Ltd. řti

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Tel/Web: +90 258 371 46 56 / www.hunkarmedikal.com

Index

1.	Product Description.....	27
2.	Product Description.....	27
3.	Device Technical Specification.....	28
4.	Device's Design and Features	29
5.	Device's Installation Instruction	30
6.	Device's Power Option.....	31
7.	Device's Usage Instruction	33
8.	Cleaning and Disinfecting	35
9.	Maintenance and Clues.....	37
10.	Device's Technical Data's	38
11.	Contraindications and Indications	39
12.	Keeping and Transportation	41
13.	Troubleshooting	42
14.	Warranty Conditions and After Sales Support.....	42
15.	EMC (Electromagnetic Compatibility) Decalaration.....	43
16.	Symbols Descriptions	48
17.	Device & Outside Box's Symbols Descriptions	49
18.	European Conformity Mark.....	49
19.	Manufacturer and Technical Service Information	49

This product is medical device. Before using the device, please read user manual carefully and use it with in the direction of your doctor.

1. Product Description

Nimo® Portable Mesh Nebulizer is designed to be used in the treatment of respiratory diseases such as asthma, COPD, chronic laryngitis, laryngitis, upper respiratory tract infections and allergies.

HNK-MESH-01 is combined with piezzo which converts electircy energy into mechanical energy and on it there are pores in micrometer size and it has active vibrating mesh technology and with the help of this technology it converts liquid medication to aerosol and provides it to take inside to body via respiration.

This device is suitable for patients from all different ages.

Active vibrating mesh technology **Nimo®** Portable Mesh Nebulizer gives the medication to patient silently, fastly and without disturbing the patient.

HNK-MESH-01 is easy to use, light and ergonomic and is designed for use at home, at work and in travelling.

2. Product Description

Active vibrating mesh technology high density aerosol exit.

Effective treatment with low MMAD value.

Low medication residual volume.

Low power consumption and silent operation.

It can be use 60 minutes with full charge.

Long life rechargeable Li-Titanate battery technology.

It can be charged in 90 minutes via USB cable.

Automatic shut-down in 10 minutes.

Easily removable medication cup design.


Light and portable mini compact design.

3. Device Technical Specification

Product Name	Portable Mesh Nebulizer
Model Number	HNK-MESH-01
Method of Operation	Active vibrating mesh technology
Power Supply	DC 5.0V 1A or Integrated Battery
Power Consumption	4.0 W
Integrated Battery	Lithium Titanate Battery 2.4 VDC 1000mAh
Nebulization Rate	≥0.3 ml/min
MMAD	<5 μm
Working Frequency	120 kHz±10 kHz
Medication Cup Capacity	10 ml
Working Environment	Temperature: 5 °C – 40 °C Humidity: <80% R.H. Atmospheric pressure: (86-106) kPa
Storage/Delivery Environment	Temperature: -20 °C – 55 °C Humidity: <80% R.H. Atmospheric pressure: (70-106) kPa
Classification	Class II Equipment
Type	Type BF Applied Parts
Ingress Protection	IP22
Noise	<40 dB
Size	41mm(L)x45mm(W)x91mm(H)
Weight	110 gr (include the battery)

4. Device's Design and Features

a) Device's appearance and pieces


Piece Number	Piece Name	Piece Number	Piece Name
1	Medication Cup Lid	6	Device On-Off Key (Led Light Indicator)
2	Medication Cup		
3	Mesh Disc	7	USB Cable
4	Medication Cup Take-Out Key	8	Adult and Child Masks
5	Device Main Unit		

b) Key Specifications

Key	Function
	On-Off
PRESS	Medication Cup Take-Out Key

c) Led Light Indicator and Meanings

Blue Light: Device is Working

Green Light: Device is Charging


No Light: Device is not Working

5. Device's Installation Instruction

Before installation, please clean and disinfect all the pieces.

- 1) First open the whole package then take-out the device and accessories.
- 2) How to set medication cup to the main unit.


As it is shown in the Picture put the medication cup in to the housing in the main unit and push down.

When you hear the **click** be sure the medication cup set correctly.

3) How to install accessories.


As it is shown in the Picture install the mask to the device's connector part and be sure it fits correctly.

6. Device's Power Option

AC/DC Charge Adaptor (Not include in the package)

- Device is without adaptor. Please use adaptor which has GB 9706.1-2007 or IEC60601-1 standards.
- Technical parameters of power adaptor should be suitable with DC.50V/1.0 A.
- Device's power system is rechargeable titanate battery.

How To Charge?

As it is shown in the picture plug the USB cable to the micro-USB connector. While the device is charging it will not continue nominal working.

- When the low battery level is first detected by the device, the blue light flashes 5 times to stop atomization, and then the device switches itself off automatically. Till the device recharge again in each turn-on the blue light flashes 5 times without atomization and it shutdown itself automatically.
- In this situation, please charge the device. Charge duration is 90 minutes in 5V/1A current.
- Use charge adaptor to charge the device. While charging the green light flashes slowly.
- When the device full charged the green light will be stable.
- With full charge the device can be used for 60 minutes.
- The device's battery life is 5 year and in a year it can charge and discharge for 3.000 cycles.
- Even you don't use the device it is advised to charge and discharge one time in a month.

Don't fix or remove the device's equipments.

Don't remove or change internal rechargeable battery. If you need to replace, ask Hünkar Eczacı Medikal.

When you finish charging the device unplug the charge adaptor.


The device and accessories which date of expiry finish should parry according to recycling and environmental protecting regulations.

7. Device's Usage Instruction

Before using the device clean and disinfect all the parts.

1) Filling the medication

Please turn
the lid
reverse side
of clockwise

Take-out the
lid by
pulling up

Put the
medication
in to cup

Close it by
turning the
lid
clockwise

Follow your doctor's advise while putting the medication and don't pass the maximum capacity of medication cup which is 10 ml.


After putting the medication to the cup in order not to pour close the lid till end of the cure.

2) Treatment

- Push the on/off  button for a second and turn-on the device and be sure the blue light is on.

- b. In the treatment process hold the device upright and be sure aerosol is coming properly from mesh.
- c. According to the user's personal needs, the user can use this device in the following way.

*** Treatment with mask**

- d. **Nimo®** Portable Mesh Nebulizer's nominal working duration is 10 minutes. After 10 minutes the device shutdowns automatically. If you want to continue to use it push the on/off  button again. Please be sure there is medication in the medication cup.
- e. After cure process finish to switch-off the device push on/off  button. Empty the remaining residue in medication cup and don't use again.

It is not recommended to use the device more than 45° oblique in the treatment.

In the nebulization process light vibration doesn't effect the usage but don't shake the device.

Follow your doctor's advice. In the cure process be calm and relax.

The medication will be deposited around the medication cup connector and mesh and it is going to effect the nebulization performance. In this situation stop the nebulization and take-out the medication cup and accessories then this residue should be cleaned according to the cleaning directions.

8. Cleaning and Disinfecting


- ⚠ **There is infection risk if medication cup doesn't clean rapidly and properly.**
- ⚠ **After every usage disinfect the medication cup, lid and masks.**
- ⚠ **Before cleaning the device if it is charging stop it and unplug USB cable.**

Cleaning and disinfection methods are down below.


Cleaning after every usage.

After every cure clean the device with warm pure water. (Under 40 °C)

- 1) Take-out the device's accessories. (Mesh, Mask and USB cable)
- 2) Empty the residual in medication cup.
- 3) Put warm (not more than 40 °C) water in to medication cup.
- 4) To clean the mesh turn-on the device for 2 minutes with warm and pure water inside.
- 5) Turn-off the device.
- 6) Then take-out the medication cup as it is shown in the pictures.

a) Push the PRESS button.

b) Take-out the medication cup up.

- 7) Clean the medication cup with warm pure water.
- 8) Dry the cleaned medication cup with gauze. Be sure all the electrodes dried totally.
- 9) Clean the main body of device with the gauze or wet cloth. Be sure electrodes are totally dried.
- 10) After cleaned all the parts be sure all the parts dried totally.

In order not to damage the mesh disc in cleaning do not use any sharp or firm things and never touch them to the mesh disc.

Don't contact the main unit with any liquid directly since this device is totally electronic device.

Don't wash the parts in dishwasher.

Don't dry the parts in microwave.

Disinfection

It is very important to disinfect the medication cup, medication cup lid, masks after every usage.

- 1) Clean the medication cup, medication cup lid, masks with warm pure water. (Not more than 40 °C)
- 2) Put the medication cup, medication cup lid and masks in to boiled distilled warm water (less than 40 °C) and wait 30 minutes than wash all the parts with distilled warm water.
- 3) Then clean all these washed parts along 1 minute with %75 pure ethyl alcohol.
- 4) Put the medication cup, medication cup lid and masks on a towel in the clean air where children can't reach to dry.
- 5) Be sure about all these parts cleaned and dried well.

Alcohol is flammable for this reason while cleaning with alcohol stay away from fire.

9. Maintenance and Clues

- Don't exposure the mesh nebulizer and its accessories to abrasive liquids and gases.
- Don't touch the mesh with sharp or perforator devices or with your hand otherwise this can cause permanent damages and the device will be unusable.
- To prevent damage and for your safety do not put any liquid to the device on treatment process.
- Do not exposure the device medication cup and accessories directly sun light.
- Protect the device, medication cup, and accessories from harsh impulse otherwise it can cause damage in the device and it can be unusable.
- Clean the device and its parts before and after every usage.
- Do not extract, fix or any change on device otherwise this can cause electric shock, fire or leakage.
- Never use the device if it is damaged or fall in to water.


If any problem occurs while using the device first check the troubleshooting. If problem continuous apply to authorized service.

10. Device's Technical Data's

- The device is manufactured in accordance with EN ISO 27427:2019 Anaesthetic and respiratory equipment - Nebulizing systems and components
- The nebulisation system data's are given in the table in below.

MMAD (Mass Median Aerodynamic Diameter)	2.5 μm
Nebulization Rate	0.61 $\pm$ 0.02 ml/min
GSD (Geometric Standard Deviation)	2.42
Particle Size Under 5 μm	>%75

- The horizontal axis is the particle size value, the value is logarithmic distribution. The left vertical axis is the cumulative percentage of the volume, corresponding to the rising trend of the curve.

- The small particle size (MMAD) for 50% cumulative particle mass indicates that the particle sizes at the aerosol from the nebulizer outlet can reach the lower airways.
- In the Nimo Mesh nebulizer, the particles in the range of 0.5-5 μm constitute approximately 75% of the aerosol output and more.
- GSD shows aerodynamic particle size distribution and is calculated as 2.42 for Nimo Mesh Nebulizer.

11. Contraindications and Indications

- **Contraindication:** means that the treatment is unsuitable.

This device is unsuitable for Pentomidin medications.

It is not suitable for pulmonary edema patients.

It is not suitable for acute lung embolism.

- **Indications:** these are the situations in which treatment is appropriate.

It is used in the treatment of respiratory tract diseases such as asthma, COPD, chronic laryngitis, laryngitis, upper respiratory tract infections, allergies, by converting liquid drugs into aerosol mist and inhaling into the body.

- **Side effects:**

No side effects have been identified as a device. The side effects of the drugs used with the device are specifically described in the instructions for use of the products used.

- **Patient population:**

The product can be used in adults, pediatrics and infants.

The device is designed for patients from every age but babies, small children and adults who need special care should use it under the control of an adult.

- **Intended User:**

This device may be used by the patient, the caregiver, or by qualified medical professionals for home treatment, under the instructions of a legally licensed medical professional such as a doctor, nurse, or therapist.

The user should also be able to understand the general operating principles of the device and the instructions in the device user manual.

- **Warnings:**

Before use the device take-out all the protective equipments.

This device is used for breathing passage treatment under the control of doctor.
Please do not misuse it except for its purpose.

If you feel yourself bad while using the device stop using and ask your doctor for help.

Please use original parts and accessories. Warranty is invalid if you use not original parts and user's personal reasons.

Before use the device be sure all the accessories are working.

To avoid the medication residue on the face, be sure to wipe the face after removing the mask.

Each user should use her/his own medication cup and mask

Don't plug or unplug the device's adaptor with wet hand.

Keep the device in the storage and transportation environment conditions defined on the user manual and product box, use it in environments with working ambient conditions.

Because of the choke hazard put the USB cable away from children.

While using the device be sure the liquid in the medication cup is touching totally to the mesh disc.

While using the device, do not keep it tilted more than 45 degrees, avoid excessive movements during treatment. Do not shake the device, protect it from impacts and falls.

Using and buying of this device should be advised by your doctor. Please follow your doctor order in medication type, limit and using way.

Do not use devices which fall in to water and don't charge it. In this situation ask help from authorised service.

12. Keeping and Transportation

- The device can be kept with its package in well aired room humid is not more than %80 and the temperature between -20 ~ 55 °C.
- Under this storage condition explained up this device can be kept for 5 years.
- While carrying protect the package from get wet, impact and being crushed.

Keep the device away from children.

Keep the device away from direct sunlight and abrasive gases.

Package Contents:

- One **Nimo®** HNK-MESH-01 portable mesh nebulizer
- User manual
- Warranty

Accessories inside the case;

- USB cable (1 pc)
- Adult mask (1 pc)
- Child mask (1 pc)

13. Troubleshooting

Problem	Solution
Device is active but not working	☞ Check the charge level. ☞ Check the medication lid sets correctly
Low aerosol or there is no aerosol	☞ Check the medication cup is filled with water soluble and true medication. ☞ Check the medication cup if there is right size of medication inside. ☞ Mesh can be clogged. Put 3-6 ml of water and 2-3 drop vinegar in to medication cup and clean the medication cup in order to use next time and disinfect it. ☞ Check the charge level. ☞ Check the medication cup sets correctly.
The blue light flashes 5 times and the device shutdown itself.	☞ The device's charge is nearly finishing. Please follow the charging instructions.
The blue light flashes 3 times and the device shutdown itself.	☞ Check the medication cup is in right place. ☞ Check there is medication in medication cup.

If these suggestions doesn't solve the problem please don't try to fix it and ask help from authorized service.

14. Warranty Conditions and After Sales Support

The device's warranty duration is for 2 years after you buy it. The warranty period of the medicine container is 6 months from the date of purchase.

The expecting life for this device is 5 years if you follow the instructions in manual charge and discharge properly.

- ! **Warning:** The warranty will be invalid if the device is open by user or harmed or try to fix it.
- ! **Note:** Please do not forget to take the warranty and sales receipts and keep them.

15. EMC (Electromagnetic Compatibility) Declaration

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual.

This device has been tested for Medical Products in accordance with EN-60601-1-2 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual.

This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment.

See the following tables for more information about this device and EMC.

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment – guidance
Emissions RF CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharges IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 0,5 kV, ± 1 kV	± 0,5 kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	± 0,5 kV, ± 1 kV, ± 2 kV	
Voltage dips IEC 61000-4-11	0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment.
	0% U_T ; 1 cycle and 70% U_T ;25/30 cycle Single phase: 0°	0% U_T ;1 cycle and 70% U_T ;25/30 cycle Single phase: 0°	
Voltage interruptions IEC 61000-4-11	0% U_T ;250/300 cycle	0% U_T ;250/300 cycle	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Not: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz 10V/m 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz 10V/m 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distance: $d = 0.35\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 1.2\sqrt{P} \quad 800 \text{ MHz to } 2.7 \text{ GHz}$ $d = 2.3\sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
f. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. g. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. h. Calibration for current injection clamps shall be performed in a 150 Ω system. i. If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range. j. Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.			

Recommended separation distances between portable and mobile RF communications equipment and the device				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands)	150 kHz to 80 MHz (out ISM and amateur radio bands)	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1.2\sqrt{P}$	$d = 0.6\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.06	0.12	0.23
0.1	0.38	0.19	0.38	0.73
1	1.2	0.6	1.2	2.3
10	3.8	1.9	3.8	7.3
100	12	6	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz. Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Guidance and manufacturer's declaration – electromagnetic emissions							
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	380-390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13.17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5240						
	5785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.
a. For some services, only the uplink frequencies are included. b. The carrier shall be modulated using a 50 % duty cycle square wave signal. c. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.
The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = 6/d\sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

16. Symbols Descriptions

	This device complies with the requirements of the Medical Devices Directive(93/42/EEC)
	Manufacturer
	Date of Manufacture
	Class II Equipment per IEC 60601-1
	Type BF Applied Part per IEC 60601-1
	Consult instructions for use
	LOT Number
IP22	Protected against foreign objects equal to or greater than 12.5mm in diameter and against drops of water falling at up to 15 °C from vertical
	Model Number
	Serial Number
	Medical Device
	Unique Device Identification
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.

17. Device & Outside Box's Symbols Descriptions

	Warning/Caution/Note
	Keep Away From Sunlight
	Keep Away From Rain
	Handle With Care
	Operating Temperature Limits: 5°C-40°C
	Storage And Delivery Temperature Limits: -20°C-55°C
	Operating, Storage And Delivery Humidity Limits: <80%R.H.
	Barcode Number

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The product conforms to European Medical Device Directive (93/42/EEC) and meets applicable health, safety and environmental requirements. If the mark is accompanied by a number, conformity is verified by the indicated notified body.

19. Manufacturer and Technical Service Information

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Unvanı : Hünkar Ecza Ve Medikal Sanayi Ve Ticaret Ltd. Şti