HNK-MESH-01

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Index

1.	Product Description.....	2
2.	Product Description.....	2
3.	Device Technical Specification.....	3
4.	Device's Design and Features	4
5.	Device's Installation Instruction	5
6.	Device's Power Option.....	6
7.	Device's Usage Instruction	8
8.	Cleaning and Disinfecting	10
9.	Maintenance and Clues.....	12
10.	Device's Technical Datas.....	13
11.	Contraindications and Warnings.....	14
12.	Keeping and Transportation	15
13.	Troubleshooting	16
14.	Warranty Conditions and After Sales Support.....	16
15.	EMC (Electromagnetic Compatibility) Declaration.....	17
16.	Symbols Descriptions	21
17.	Device & Outside Box's Symbols Descriptions	22
18.	European Conformity Mark.....	23
19.	Manufacturer and Technical Service Information	23

This product is medical device. Before using the device please read user manual carefully and use it with in the direction of your doctor.

1. Product Description

Nimo® Portable Mesh Nebulizer is designed for asthma, allergy and other respiratory disease.

HNK-MESH-01 is combined with piezoelectric which converts electrical energy into mechanical energy and on it there are pores in micrometer size and it has titanium alloy active mesh technology and with the help of this technology it converts liquid medication to aerosol and provides it to take inside to body via respiration.

This device is suitable for patients from all different ages.

Titanium alloy active mesh technology NIMO Portable Mesh Nebulizer gives the medication to patient silently, fastly and also without disturbing the patient.

HNK-MESH-01 is easy to use, light and ergonomic and is designed for use at home, at work and in travelling.

2. Product Description


- With titanium alloy mesh technology high density aerosol exit.
- 2.5 µm MMAD efficient treatment.
- Low medication residual volume.
- Low power consumption and silent operation.
- It can be used 60 minutes with full charge.
- Long life rechargeable Li-Titanate battery technology.
- It can be charged in 90 minutes via USB cable.
- Automatic shut-down in 10 minutes.
- Automatic shut-down when medication finishes.
- Easy removable medication cup design.
- Light and portable mini compact design.

3. Device Technical Specification

Product Name	Portable Mesh Nebulizer
Model Number	HNK-MESH-01
Method of Operation	Active vibrating mesh technology
Power Supply	DC 4.8V (Inbuilt lithium titanate battery) or DC 5.0V 1A
Power Consumption	<4.0 W
Nebulization Rate	0.5±0.15 ml/min
MMAD	2.5±30%µm
Working Frequency	120 kHz±10 kHz
Medication Cup Capacity	10 ml
Working Environment	Temperature: 5 °C – 40 °C Humidity: <80% R.H. Atmospheric pressure: (86.0-106.0)kPa
Storage/Delivery Environment	Temperature: -20 °C – 55 °C Humidity: <80% R.H. Atmospheric pressure: (70-106.0)kPa
Classification	Class II Equipment
Type	Type BF Applied Parts
Ingress Protection	IP22
Noise	<5 dB
Size	41mm(L)x45mm(W)x91mm(H)
Weight	110 gr (include the battery)

4. Device's Design and Features

a) Device's appearance and pieces


Piece Number	Piece Name	Piece Number	Piece Name
1	Medication Cup Lid	6	Device On-Off Key (Led Light Indicator)
2	Medication Cup		
3	Titanium Alloy Mesh	7	USB Cable
4	Medication Cup Take-Out Key	8	Adult and Child Masks
5	Device Main Unit	9	Mouth Piece

b) Key Specifications

Key	Function
	On-Off
PRESS	Medication Cup Take-Out Key

c) Led Light Indicator and Meanings


Blue Light: Device is Working
Green Light: Device is Charging
No Light: Device is not Working

5. Device's Installisation Instruction

Before installisation please clean and disinfection all the pieces.

- 1) First open the whole package then take-out the device and accessories.
- 2) How to set medication cup to the main unit.


As it is shown in the Picture put the medication cup in to the nest in the main unit and push down.


When you hear the **click** be sure the medication cup set correctly.

3) How to install accessories.

As it is shown in the Picture install the mask to the device's connector part and be sure it fits correctly.


As it is shown in the Picture install the mouth piece to the device's connector part and be sure it fits correctly.

6. Device's Power Option

AC/DC Charge Adaptor (Not include in the package)

- Device is without adaptor. Please use adaptor which has GB 9706.1-2007 or IEC60601-1: 2012 standarts
- Technical parameters of power adaptor should be suitable with DC.50V/1.0 A.
- Device's power system is rechargeable titanet battery.

How To Charge?

As it is shown in the picture plug the USB cable to the mini USB connector. While the device is charging it continues nominal working.

- When the low battery level is first detected by the device, the blue light flashes 5 times to stop atomization, and then the device switches itself off automatically. Till the device recharge again in each turn-on the blue light flashes 5 times without atomisation and it shut-down itself automatically.
- In this situation please charge the device. Charge duration is 90 minutes in 5V/1A current.
- Use charge adaptor to charge the device. While charging the green light flashes slowly.
- When the device full charged the green light will be stable.
- With full charge the device can be used for 60 minutes.
- The device's battery life is 5 year and in a year it can charge and discharge for 3.000 cycles.
- Even you dont use the device it is advised to charge and discharge one time in a month.

Dont fix or remove the device's equipments.

Dont remove or change internal rechargeable battery. If you need to chance ask Hünkar Medikal ve Eczacı.


When you finish the charging the device unplug the charge adaptor.

 The device and accessories which date of expiry finish should comply according to recycling and environment protecting regulations.

7. Device's Usage Instruction

 Before using the device clean and disinfect all of the parts.

1) Fullfill Medication

Please turn
the lid
reverse side
of clockwise

Take-out the
lid by
pulling up

Put the
medication
in to cup

Close it by
turning the
lid
clockwise


 Follow your doctor's advice while putting the medication and don't pass the maximum capacity of medication cup which is 10 ml.

 After putting the medication to the cup in order not to pour close the lid till end of the cure.

2) Cure


- a. Push the on/off  button for a second and turn-on the device and be sure the blue light is on..
- b. In the cure process hold the device upright and be sure aerosol is coming properly from mesh.
- c. According to user's personal needs user can use this device in these two different ways.

*Cure with mask

As it is shown in the picture put the mask on to cover mouth and nose and breathe normally.
Device's Usage Instruction

* Cure with mouth piece

Put the mouth piece as it is shown in the picture between your lips and tooth and breathe from your mouth normally.

- d. **Nimo**® Portable Mesh Nebulizer's nominal working duration is 10 minutes. After 10 minutes the device shut-downs automatically. If you want to continue to use it push the on/off  button again. Please be sure there is medication in the medication cup.
- e. After cure process finish to switch-off the device push on/off  button. Empty the

remaining residue in medication cup and don't use again.

It is not recommended to use the device more than 45° oblique in the cure.

In the nebulization process light vibration doesn't affect the usage but don't shake the device.

Follow your doctor's advice. In the cure process be calm and relax.

The medication will be deposited around the medication cup connector and mesh and it is going to affect the nebulization performance. In this situation stop the nebulization and take-out the medication cup and accessories then this residue should be cleaned according to the cleaning directions.

8. Cleaning and Disinfecting

There is infection risk if medication cup doesn't clean rapidly and properly.

After every usage disinfect the medication cup, lid, mouthpiece and masks.

Before clean the device if it is charging stop it and unplug USB cable..


Cleaning and disinfection methods are down side.

Cleaning after every usage.


After every cure clean the device with warm pure water. (Under 40 °C)

- 1) Take-out the device's accessories. (Mesh, mouth piece, usb cable)
- 2) Empty the residual in medication cup.
- 3) Put warm (not more than 40 °C) water in to medication cup.
- 4) To clean the mesh turn-on the device for 2 minutes with warm and pure water inside.
- 5) Turn-off the device.

6) Then take-out the medication cup as it is shown in the pictures.

a) Push the PRESS button.

b) Take-out the medication cup up.

7) Clean the medication cup with warm pure water.

8) Dry the cleaned medication cup with gauze. Be sure all the electrodes dried totally.

9) Clean the main body of device with the gouze or wet cloth. Be sure electrodes are totally dried.

10) After cleaned all the parts be sure all the parts dried totally.

In order not to damage the mesh disc in cleaning do not use any sharp or firm things and never touch them to the mesh disc.

Don't contact the main unit with any liquid directly since this device is totally electronic device.

Don't wash the parts in dishwasher.

Don't dry the parts in microwave.

Disinfection

It is very important to disinfect the medication cup, medication cup lid, mouth piece, masks after every usage.

- 1) Clean the medication cup, medication cup lid, mouth piece, masks with warm pure water. (not more than 40 °C)
- 2) Put the medication cup, medication cup lid, mouth piece and masks in to boiled distilled warm water (less than 40 °C) and wait 30 minutes then wash all the parts with distilled warm water.
- 3) Then clean all these washed parts along 1 minute with %75 pure ethyl alcohol.
- 4) Put the medication cup, medication cup lid, mouth piece and masks on a towel in the clean air where children can't reach to dry.
- 5) Be sure about all these parts cleaned and dried well.

Alcohol is flammable for this reason while cleaning with alcohol stay away from fire.

9. Maintenance and Cues

- Don't expose the mesh nebulizer and its accessories to abrasive liquids and gases.
- Don't touch the mesh with sharp or perforator devices or with your hand otherwise this can cause permanent damages and the device will be unusable.
- In order to prevent damage and for your safety do not put any liquid to the device during cure process.
- Do not expose the device medication cup and accessories directly to sunlight.
- Protect the device, medication cup, and accessories from harsh impulse otherwise it can cause damage in the device and it can be unusable.
- Clean the device and its parts before and after every usage.
- Do not extract, fix or any change on device otherwise this can cause electric shock, fire or leakage.
- Never use the device if it is damaged or fall in to water.


If any problem occurs while using the device first check the troubleshooting. If problem continues apply to authorized service.

10. Device's Technical Datas

- The device is manufactured in compliance with EN 13544-1: 2007 + A1: 2009 standards for airway treatment devices.
- The nebulisation system datas are given in the table in down side.

MMAD (Mass Median Aerodynamic Diameter)	2.5±30%µm
Nebulization Rate	0.5 ± 0.15 ml/min
GSD (Geometric Standard Deviation)	2.42
Particle Size Under 5 µm	>%75

- The horizontal axis is the particle size value, the value is logarithmic distribution; The left vertical axis is the cumulative percentage of the volume, corresponding to the rising trend of the curve.

- The small particle size (MMAD) for 50% cumulative particle mass indicates that the particle sizes at the aerosol from the nebulizer outlet can reach the lower airways.

- In the Nimo Mesh nebulizer, the particles in the range of 0.5-5 μm constitute approximately 75% of the aerosol output and more.
- GSD shows aerodynamic particle size distribution and is calculated as 2.42 for Nimo Mesh Nebulizer.

11. Contraindications and Warnings

* **Contraindication means that the cure is unsuitable.**

- This device is unsuitable for Pentomidin medications.
- It is not suitable for pulmonary edema patients.
- It is not suitable for acute lung embolism.

* **Warnings**

The device is designed for patients from every age but babies, small children and adults who need special care should use it under the control of an adult.

Before use the device take-out all the protective equipments.

This device is used for breathing passage treatment under the control of doctor. Please do not misuse it except for its purpose.

If you feel yourself bad while using the device stop using and ask your doctor for help.

Please use original parts and accessories. Warranty is invalid if you use not original parts and user's personal reasons.

Before use the device be sure all the accessories are working.

To avoid the medication residue on the face, be sure to wipe the face after removing the mask.

Each user should use her/his own medication cup, mask and mouth piece.

Don't plug or unplug the device's adaptor with wet hand.

Because of the drowning danger put the USB cable away from children.

While using the device be sure the liquid in the medication cup is touching totally to the mesh disc.

Using and buying of this device should be advised by your doctor. Please follow your doctor order in medication type, limit and using way.

Do not use devices which fall in to water and don't charge it. In this situation ask help from authorised service.

12. Keeping and Transportation

- The device can be kept with its package in well aired room humid is not more than %80 and the tempeture between -25~50 °C.
- Under this storage condition explained up this device can be kept for 5 years.
- While carrying protect the package from get wet, impact and being crushed.

Keep the device away from children.

Keep the device away from direct sunlight and abrasive gases.

Package Contents

- One HNK-MESH-01 Nebulizer.
- User's Manual
- Warranty

Accessories inside the case

- USB Cable (1)
- Adult Mask (1)
- Child Mask (1)
- Mouth Piece (1)

13. Troubleshooting

Problem	Solution
Device is active but not working	<ul style="list-style-type: none"> ☞ Check the charge level. ☞ Check the medication lid sets correctly
Low aerosol or there is no aerosol	<ul style="list-style-type: none"> ☞ Check the medication cup is filled with water soluble and true medication. ☞ Check the medication cup if there is right size of medication inside. ☞ Mesh can be cöngest. Put 3-6 ml of water and 2-3 drop vinegar in to medication cup and clean the medicaiton cup in order to use next time and disinfection it. ☞ Check the charge level. ☞ Check the medication cup sets correctly.
The blue light flashes 5 times and the device shut-down itself.	<ul style="list-style-type: none"> ☞ The device's charge is nearly finishing. Please follow the charging instructions.
The blue light flashes 3 times and the device shut-down itself.	<ul style="list-style-type: none"> ☞ Check the medication cup is in right place. ☞ Check there is medication in medication cup.

If these suggestions doesn't solve the problem please don't try to fix it and ask help form authorized service.

14. Warranty Contitions and After Sales Support

The device's warranty duration is for 2 years after you buy it. The warranty period of the medicine container is 6 months from the date of purchase.

The expecting life fort his device is 5 years if you follow the instructions in manual charge and decharge properly.

- ! **Warning: The warranty will be invalid if the device is open by user or harmed or try to fix it.**
- ! **Note: Please do not forget to take the warranty and sales receipts and keep them.**

15. EMC (Electromagnetic Compatibility) Declaration

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual.

This device has been tested for Medical Products in accordance with EN-60601-1-2 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual.

This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment.

See the following tables for more information about this device and EMC.

Guidance and Manufactures Declaration-Electromagnetic		
These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment—guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufactures Declaration-Electromagnetic			
These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.			
Immunity Test	IEC60601 Test Level	Compliance Level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	± 6 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential Mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<%5 UT (UT ^{de} >%95 düşüş) 0,5 devir için %40 UT (UT ^{de} %60 düşüş) 5 devir için %70 UT (UT ^{de} %30 düşüş) 25 devir için <%5 UT (UT ^{de} >%95 düşüş) 5 saniye için	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power source
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: Ut is the AC mains voltage prior to application of the test level.			

Guidance and Manufactures Declaration-Electromagnetic

These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.

Immunity Test	IEC60601 Test Level	Compliance Level	Electromagnetic environment—guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V rms 150 kHz ila 80 MHz 10 V rms 150 kHz ila 80 MHz 10 V/m 80 MHz ila 2,5 GHz	3 V 10 V 10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = 1,16\sqrt{P}$ $d = 1,20\sqrt{P}$ $d = 1,2\sqrt{P} \text{ 80 MHz ila 800 MHz}$ $d = 2,3\sqrt{P} \text{ 800 MHz ila 2,5 GHz}$ <p>Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Not 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

d. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio,AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

e. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distance between portable and mobile RF communication equipment and this device

These devices are intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz ila 80 MHz, Outside ISM Bands $d = 1,16\sqrt{P}$	150 kHz ila 80 MHz, In ISM Bands $d = 1,20\sqrt{P}$	80 MHz ila 800 MHz $d = 4\sqrt{P}$	80 MHz ila 800 MHz $d = 7,66\sqrt{P}$
0,01	0,12	0,12	0,12	0,23
0,1	0,37	0,38	0,38	0,73
1	1,16	1,20	1,20	2,30
10	3,67	3,79	3,79	7,27
100	11,60	12,00	12,00	23,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies..



Note 2. ISM (industrial, scientific, and medical) bands from 150 kHz to 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz.

Note 3. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

16. Symbols Descriptions

	This device complies with the requirements of the Medical Devices Directive(93/42/EEC)
	Manufacturer
	Date of Manufacture
	Class II Equipment per IEC 60601-1
	Type BF Applied Part per IEC 60601-1
	Consult instructions for use
	LOT Number
	Protected against foreign objects equal to or greater than 12.5mm in diameter and against drops of water falling at up to 15 °C from vertical
	Model Number
	Serial Number
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.

17. Device & Outside Box's Symbols Descriptions

	Warning/Caution/Note
	Keep Away From Sunlight
	Keep Away From Rain
	Handle With Care
	Operating Temperature Limits: 5°C-40°C
	Storage And Delivery Temperature Limits: -20°C-55°C
	Operating, Storage And Delivery Humidity Limits: <80%R.H.
	RoHS Declaration
	Reach Declaration
	Barcode Number

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The product conforms to European Medical Device Directive (93/42/EEC) and meets applicable health, safety and environmental requirements. If the mark is accompanied by a number, conformity is verified by the indicated notified body.

19. Manufacturer and Technical Service Information

Hünkar Ecza ve Medikal San. Tic. Ltd. Şti.
Akçeşme Mahallesi 2019 sk. No:5 Merkezefendi
Denizli /TÜRKİYE
Tel: +90 258 371 46 56
www.hunkarmedikal.com