

nimo®

STETESKOP



Yüksek kalite, Daha iyi sağlık hizmeti

TEK BAŞLI STETESKOP

Steteskop insan vücudundaki kalp,ciğer ve diğer sesleri duymak için üretilmiştir.

Steteskop; hafif dayanaklı aliminyum diyafram, ayarlanabilir krom kaplı binaural,“Y” şeklinde PVC boru ve beyaz ABS kulak tıkaçları ile birlikte üretilmiştir. Ekseninde dönebilen başlığı da mevcuttur.Ayrıca steteskopun kullanımı kolay,kaliteli ve ses kaliteside oldukça yüksektir.



Manufacturer / Üretici Firma: Honsun (Nantong) CO.Ltd.
No:8 Tongxing Road Economic & Techonogical
Development Area 226009 Nantong City P.R.C
www.asian-medical.com
+862163056696

Importer / İthalatçı Firma: Hünkar Ecza ve Medikal Ltd. Şti.
Akçeşme Mevki Sk. 2019 No/5 Merkezefendi/ DENİZLİ-TURKEY
+90 258 371 46 56
www.hunkarmedikal.com

nimo®

STETHOSCOPE



High quality, Better healthcare

SINGLE HEAD STETHOSCOPE

Used for listening various sounds of heat,lung and organs in human body.

The single head model consists of lightweight aluminum chest-piece, adjustable chrome-plated brass binaural, durable one-piece molded PVC “Y” tubing and pure white ABS eartips. Turning the convertible axis to function each head.models provide good acoustic, reponse and comfort fit.

